

S. aureus Gyrase Microplate Assay Kit



Product Description (Product Number # SATRG01, SATRG02)

The kit is supplied with sufficient *S. aureus* gyrase enzyme, plasmid DNA substrate (relaxed pNO1; supplied at 1 mg/ml), 5X Assay Buffer, Enzyme Dilution Buffer and TFO1 oligo for 100 assays. The enzyme is supplied at a concentration of 10 U/ μ l in Dilution Buffer. The kit is also supplied with sufficient Wash Buffer, TF buffer and T10 buffer for one 96-well plate. These buffers are supplied as 20X concentrates and must be diluted prior to use with ultra-pure water.

Store at -80 °C. (Stable for 3 months undiluted.) It is recommended that the enzyme is aliquoted to avoid repeated freeze-thaw cycles.

For *in vitro* laboratory research use only.

Dilution Buffer

50 mM Tris.HCl (pH 7.5)
1 mM DTT
1 mM EDTA
40 % (w/v) glycerol

Assay Buffer (supplied as 5X stock)

40 mM HEPES.KOH (pH 7.6)
10 mM magnesium acetate
10 mM DTT
2 mM ATP
500 mM potassium glutamate
0.05 mg/ml albumin

TF Buffer (supplied as a 20X stock)

50 mM sodium acetate (pH 5.0)
50 mM NaCl
50 mM MgCl₂

Wash Buffer (supplied as a 20X stock)

20 mM Tris.HCl (pH 7.6)
137 mM NaCl
0.01 % (w/v) BSA
0.05 % (v/v) Tween-20

T10 Buffer (supplied as a 20X stock)

10 mM Tris.HCl (pH 8)
1 mM EDTA

Preparation of Plate and Supercoiling Assay

Rehydrate wells with 3 x 200 μ l Wash Buffer (diluted from 20X stock before use).

Immobilize 100 μ l of 500nM TFO1 oligo in each well (5 μ l of 10 μ M TFO1 in 95 μ l Wash Buffer), 5 minutes at room temperature. Wash off excess oligo with 3 x 200 μ l Wash Buffer.


Incubate 1.5 U of *S. aureus* gyrase with 0.75 μ g of relaxed pNO1 in a reaction volume of 30 μ l at 37°C for 30 minutes in Assay Buffer. Incubate reaction in well.

Add 100 μ l TF Buffer (diluted from 20X stock before use) to well and incubate for a further 30 minutes at room temperature to allow triplex formation.

Remove liquid from well and wash with 3 x 200 μ l TF Buffer to remove unbound plasmid.

Stain with appropriate fluorescence stain (Suggested stain, SYBR Gold[®] (Invitrogen) diluted to 1X with T10 buffer. Add 200 μ l per well. Incubate for 10 - 20 minutes, mix and read in fluorescence plate reader; Ex: 495 nm; Em: 537 nM).

Supercoiling of relaxed pNO1 by varying amounts of *S. aureus* gyrase in the plate assay



Quality Control

- 1) Purity: *S. aureus* gyrase is purified to > 95 % purity as judged by SDS-polyacrylamide gel electrophoresis.
- 2) pBR322 was also incubated for 4 hrs in assay buffer + 10 mM MgCl₂ at 37 °C. These tests were negative for the formation of linear products, indicating the absence of nuclease contamination.
- (3) No activity was detectable when the single subunits were assayed alone.

References

Maxwell, A., Burton, N.P. and O'Hagan, N. (2006) High-throughput assays for DNA gyrase and other topoisomerases. *Nucleic Acid Res.* **34(15)**, e104

Patent held by Plant Biosciences Ltd., Norwich, Norfolk, UK (Patent No. GB0424953.8.)

Kit issued with limited licence for individual use only