

S. aureus Topo IV Microplate Assay Kit



Product Description (Product Number # SATRIV01, SATRIV02)

The kit is supplied with sufficient *S. aureus* topoisomerase IV enzyme, plasmid DNA substrate (supercoiled pNO1; supplied at 1 mg/ml), 5X Assay Buffer, Enzyme Dilution Buffer and TFO1 oligo for 100 assays. The enzyme is supplied at a concentration of 10 U/μl in Dilution Buffer. The kit is also supplied with sufficient Wash Buffer, TF buffer and T10 buffer for one 96-well plate. These buffers are supplied as 20X concentrates and must be diluted prior to use with ultra-pure water.

Store at -80 °C. (Stable for 3 months undiluted.) It is recommended that the enzyme is aliquoted to avoid repeated freeze-thaw cycles.

For *in vitro* laboratory research use only.

Dilution Buffer

50 mM Tris.HCl (pH 7.5)
1 mM DTT
1 mM EDTA
40 % (w/v) glycerol

Assay Buffer (supplied as 5X stock)

50 mM Tris.HCl (7.5)
5 mM MgCl₂
5 mM DTT
1.5 mM ATP
350 mM potassium glutamate
0.05 mg/ml albumin

TF Buffer (supplied as a 20X stock)

50 mM sodium acetate (pH 5.0)
50 mM NaCl
50 mM MgCl₂

Wash Buffer (supplied as a 20X stock)

20 mM Tris.HCl (pH 7.6)
137 mM NaCl
0.01 % (w/v) BSA
0.05 % (v/v) Tween-20

T10 Buffer (supplied as a 20X stock)

10 mM Tris-HCl (pH 8)
1 mM EDTA

Preparation of Plate and Relaxation Assay

Rehydrate wells with 3 x 200 μl Wash Buffer (diluted from 20X stock before use).

Immobilize 100 μl of 500nM TFO1 oligo in each well (5 μl of 10 μM TFO1 in 95 μl Wash Buffer), 5 minutes at room temperature. Wash off excess oligo with 3 x 200 μl Wash Buffer.


Incubate 1.5 U of *S. aureus* topo IV with 0.75 μg of supercoiled pNO1 in a reaction volume of 30 μl at 37°C for 30 minutes in Assay Buffer. Incubate reaction in well.

Add 100 μl TF Buffer (diluted from 20X stock before use) to well and incubate for a further 30 minutes at room temperature to allow triplex formation.

Remove liquid from well and wash with 3 x 200 μl TF Buffer to remove unbound plasmid.

Stain with appropriate fluorescence stain (Suggested stain, SYBR Gold[®] (Invitrogen) diluted to 1X with T10 buffer. Add 200 μl per well. Incubate for 10 - 20 minutes, mix and read in fluorescence plate reader; Ex: 495 nm; Em: 537 nM).

Relaxation of supercoiled pNO1 by varying amountsof *S. aureus* Topo IV in the plate assay



Quality Control

1) Purity: *S. aureus* topo IV is purified to > 95 % purity as judged by SDS-polyacrylamide gel electrophoresis. 2) pBR322 was also incubated for 4 hrs in assay buffer + 10 mM MgCl₂ at 37 °C. These tests were negative for the formation of linear products, indicating the absence of nuclease contamination. (3) No activity was detectable when the single subunits were assayed alone.

References

Maxwell, A., Burton, N.P. and O'Hagan, N. (2006) High-throughput assays for DNA gyrase and other topoisomerases. *Nucleic Acid Res.* **34(15)**, e104

Patent held by Plant Biosciences Ltd., Norwich, Norfolk, UK (Patent No. GB0424953.8.)

Kit issued with limited licence for individual use only