

# *E. coli* DNA Gyrase Cleavage Assay Kit



## Product Description

(Product Numbers GCK001, GCK002, GCK003 and GCK004)

*E. coli* DNA gyrase is prepared from the overproducing strains JMtacA and JMtacB (Hallett *et al.*, 1990) and is supplied as an A<sub>2</sub>B<sub>2</sub> complex. The enzyme is supplied at a concentration of 3.5 µM in Dilution Buffer and is suitable for cleavage assays. Cleavage activity is a minimum of 2.5 – 5.0 U/µl. 50 % cleavage can be obtained with 0.2 µl in the presence of 10 µM CFX in a 30 µl reaction (see typical titration below).

Store at -80°C. For *in vitro* laboratory research use only.

### Dilution Buffer


50 mM Tris.HCl (pH 7.5)  
100 mM KCl  
2 mM DTT  
1 mM EDTA  
50 % (w/v) glycerol

### Assay Buffer (supplied as 5X stock)

35 mM Tris.HCl (pH 7.5)  
24 mM KCl  
4 mM MgCl<sub>2</sub>  
2 mM DTT  
1.8 mM spermidine  
6.5 % (w/v) glycerol  
0.1 mg/ml albumin

## Cleavage Assay

DNA gyrase is incubated with 0.5 µg of supercoiled pBR322 in a reaction volume of 30 µl at 37 °C for 1 hour in Assay Buffer in the presence of CFX. 0.2 % SDS and 0.1 mg/ml Proteinase K are added before a further incubation at 37 °C for 30 minutes.



## Quality Control

Purity: The A and B subunits are purified to > 95 % purity as judged by SDS-polyacrylamide gel electrophoresis.

Endonuclease assay: 0.5 µg relaxed pBR322 incubated with 0.5 U of DNA gyrase for 1 hour at 37 °C in the presence of 1 mM ATP shows no detectable conversion of supercoiled DNA to either open circular or linear forms when assayed by agarose gel electrophoresis.

## Reference

Hallett, P., Grimshaw, A.J., Wigley, D.B. and Maxwell, A. (1990). Cloning of the DNA gyrase genes under *tac* promoter control: overproduction of the gyrase A and B proteins. *Gene* 93: 139-142