Topoisomerase IV Cleavage Kit

(isolated from Escherichia coli)





Product Description (Product Numbers #TCK001 and TCK002)

Topoisomerase IV is prepared by overexpressing the par C and par E subunits in *E. coli* and purifying them by methods adapted from Peng and Marians, 1999. It is supplied as a heterotetramer complex. The enzyme is supplied at a concentration of 10 U/ μ I in Dilution Buffer. (6.0 μ M) 50 % cleavage can be obtained with 0.1 μ I (20 nM final concentration) in the presence of 50 μ M Norfloxacin in a 30 μ I reaction (see typical titration below).

Store at -80°C. (Stable for 12 months undiluted) For in-vitro laboratory research use only.

Dilution Buffer

40 mM HEPES-KOH (pH 7.6) 100 mM Potassium Glutamate 1 mM dithiothreitol 1 mM EDTA 40% (v/v) glycerol

Assay Buffer (supplied as 5X stock)

40 mM HEPES-KOH (pH 7.6) 100 mM Potassium Glutamate 10 mM Magnesium Acetate


10 mM dithiothreitol

50 ug/ml albumin

Cleavage Assay

1U *E. coli* Topo IV is incubated with 0.5 μg of supercoiled pBR322, in the presence of norfloxacin, in a reaction volume of 30 μl at 37°C for 1 hour in Assay Buffer.

0.2% SDS and 0.1mg/ml Proteinase K are added before a further incubation at 37° C for 30 minutes. Cleavage is detected by the presence of a linear band.



Quality Control

Purity: The parC and parE subunits are purified to >95% purity as judged by SDS-polyacrylamide gel electrophoresis.

No activity was detectable when the single subunits were assayed alone

pBR322 was also incubated for 4 hrs in assay buffer + 10 mM MgCl₂ at 37 °C. These tests were negative for the formation of linear products, indicating the absence of nuclease contamination.

References

Peng, H. and Marians, K.J. (1999) Over-expression and purification of bacterial topoisomerase IV, in DNA Topoisomerase Protocols Vol. I (Bjornsti, M-A., and Osheroff, N. eds.), Humana Press, Totowa, N. Jersey p.163-169

Marians, K.J. and Hiasa, H. (1997) Mechanism of quinolone action. A drug-induced structural perturbation of the DNA precedes strand cleavage by topoisomerase IV. J. Biol. Chem., 272, 9401-9409.

Nurse, P., Bahng, S., Mossessova, E. and Marians, K.J. (2000) Mutational analysis of Escherichia coli topoisomerase IV. II. ATPase negative mutants of ParE induce hyper-DNA cleavage. J. Biol. Chem., 275, 4104-4111.