

# Human Topo II beta Decatenation Assay Kit



## Product Description (Product Number HTKB201, HTKB202)

Human topoisomerase II beta is prepared by overexpressing in baculovirus-infected insect cells (*Spodoptera frugiperda*) and purifying it by methods developed in-house. The enzyme is supplied at a minimum concentration of 10 U/ $\mu$ l in Dilution Buffer.

Store at -80 °C.

It is recommended that the enzyme is aliquoted to avoid repeated freeze-thaw cycles.

**For *in vitro* laboratory research use only.**

### Dilution Buffer

50 mM Tris.HCl (pH 7.5)  
100 mM NaCl  
1 mM DTT  
0.5 mM EDTA  
50 % (v/v) glycerol  
50  $\mu$ g/ml albumin

### Assay Buffer (supplied as 10X stock)

50 mM Tris.HCl (pH7.5)  
125 mM NaCl  
10 mM MgCl<sub>2</sub>  
5 mM DTT  
100  $\mu$ g/ml albumin

### ATP (30X stock)


30 mM ATP

## Decatenation Assay

A typical reaction will contain 3  $\mu$ l of (10x) Assay Buffer, 1  $\mu$ l of (30x) ATP, 2  $\mu$ l of kDNA (100 ng/ $\mu$ l) plus human topoisomerase II, in a total volume of 30  $\mu$ l.

1 U of human topoisomerase II will decatenate 200 ng of kDNA when incubated in 1X Assay Buffer plus 1 mM ATP in a total reaction volume of 30  $\mu$ l at 37 °C for 30 minutes. Reactions are stopped by the addition of an equal volume chloroform/isoamyl alcohol (24:1 v:v) and another volume of 2X stop dye (40% sucrose, 10 mM EDTA, 100 mM Tris.HCL (7.5), 0.5  $\mu$ g/ml bromophenol blue).

Gels can be run in the presence or absence of ethidium bromide or chloroquine (CQ).



## Quality Control

1) Purity: Human topoisomerase II is purified to > 95 % purity as judged by SDS-polyacrylamide gel electrophoresis. 2) Tests for human topoisomerase I contamination by looking for relaxation of sc pBR322 under topoisomerase I assay conditions were negative. 3) kDNA or pBR322 were also incubated for 4 hrs in assay buffer (+ 10 mM MgCl<sub>2</sub>) at 37 °C. These tests were negative for the formation of linear products, showing the absence of nuclease contamination.

## Inhibition of decatenation activity of human topoisomerase II

The inhibition of the decatenation activity of human topo II by a compound can be determined by performing the assays as described above and including a range of concentrations of the test compound. Sufficient enzyme is included to just give full decatenation in the absence of inhibitor. Initially a wide range of concentrations of inhibitor are tested and then after determining an approximation of the  $IC_{50}$ , a range of 10X to 0.1X the  $IC_{50}$  is tested. The  $IC_{50}$  for inhibition of decatenation can be visually assessed as the concentration of compound which leads to a 50% reduction in the amount of mini-circles produced. This is then verified using gel documentation software and statistical analysis. In the example below the inhibition of human topo II-catalysed decatenation of kDNA by etoposide is followed.

