

## *M. tuberculosis* Gyrase (HIS) Decatenation kit



### Product Description (Product Numbers MTD001 and MTD002.)

*M. tuberculosis* gyrase is a C-terminal HIS-tagged protein purified after heterologous expression in *E. coli* and is supplied as an A<sub>2</sub>B<sub>2</sub> complex. The enzyme is supplied at a minimum activity of 2 U/μl in Dilution Buffer. Store at -80°C. It is recommended that the enzyme is aliquoted to avoid repeated freeze-thaw cycles.

For *in vitro* laboratory research use only.

#### Dilution Buffer (1X)

50 mM Tris.HCl (pH 7.9)  
5 mM DTT  
30 % (w/v) glycerol

#### Assay Buffer (supplied as 5X stock)

40 mM Tris.HCl (7.5)  
10 mM NaCl  
10 mM DTT  
1 mM ATP  
250 mM potassium glutamate  
6 mM Magnesium Acetate  
0.5 mg/ml albumin

#### kDNA

Supplied at 100ng/μl  
in 10mM Tris.HCl  
(pH8.0), 1mM EDTA.


The kDNA should be  
stored at 4°C.

#### Decatenation Assay

1 U of *M. tuberculosis* gyrase will decatenate 200 ng of kDNA when incubated in 1X assay buffer in a total reaction volume of 30 μl at 37°C for 30 minutes.

Gels can be run in the presence or absence of Ethidium Bromide or Chloroquine which will resolve nicked OC DNA from relaxed.

Mtub gyrase decatenation results in a range of DNA forms which migrate to form 2 distinct bands of relaxed and supercoiled DNA.



#### Quality Control

1) Purity: *M. tuberculosis* DNA gyrase is purified to > 95 % purity as judged by SDS-polyacrylamide gel electrophoresis.

2) No activity was detectable when the single subunits were assayed alone.

3) kDNA is purified to give a ratio of OD260 to OD280 ratio of 1.8 -1.9 with no detectable free mini-circles when catenated kDNA is run on a 1 % agarose gel. Each batch is tested for decatenation with purified topoisomerase IV.