

S. aureus DNA Gyrase Cleavage Assay Kit



Product Description (Product Numbers SAGC001 and SAGC002)

S. aureus DNA gyrase was purified after heterologous expression in *E. coli* and is supplied as an A₂B₂ complex. The enzyme is supplied at a concentration of 4.0 μM in Dilution Buffer and is suitable for cleavage assays. Cleavage activity is 2 U/μl.

Up to 50 % cleavage can be obtained with 0.5 μl in the presence of 10 μM CFX in a typical 30 μl reaction (see typical titration below).

Store at -80°C.

For *in vitro* laboratory research use only.

Dilution Buffer

50 mM Tris.HCl (pH 7.5)
1 mM DTT
1 mM EDTA
40 % (w/v) glycerol


Assay Buffer (supplied as 5X stock)

40 mM HEPES. KOH (pH 7.6)
10 mM magnesium acetate
10 mM DTT
100 mM potassium glutamate
0.05 mg/ml albumin

Cleavage Assay

S. aureus DNA gyrase is incubated with 0.5 μg of supercoiled pBR322 in a reaction volume of 30 μl at 37 °C for 1 hour in Assay Buffer in the presence of CFX.

0.2 % SDS and 0.1 mg/ml Proteinase K are added before a further incubation at 37 °C for 30 minutes.



Quality Control

- 1) Purity: *S. aureus* DNA gyrase is purified to > 95 % purity as judged by SDS-polyacrylamide gel electrophoresis.
- 2) pBR322 was also incubated for 4 hrs in assay buffer + 10 mM MgCl₂ at 37 °C. These tests were negative for the formation of linear products, indicating the absence of nuclease contamination.
- 3) No activity was detectable when the single subunits were assayed alone.