

# Yeast Topoisomerase II (*Saccharomyces cerevisiae*)



## Product Description (Product Number YT201 and YT205)

Yeast topoisomerase II is prepared by overexpression in *Saccharomyces cerevisiae* and purified by methods developed in-house. The enzyme is supplied at a minimum concentration of 10 U/μl in Dilution Buffer.

Store at -80 °C. (Stable for 3 months undiluted.)

**For *in vitro* laboratory research use only.**

### Dilution Buffer

50 mM Tris.HCl (pH 7.5)  
200 mM NaCl  
5 mM DTT  
1 mM EDTA  
50 % (v/v) glycerol

### Assay Buffer (supplied as 5X stock)


20 mM Tris.HCl (pH7.9)  
200 mM KCl  
10 mM MgCl<sub>2</sub>  
4 % (v/v) glycerol

### ATP (30X stock)

30 mM ATP

## Relaxation Assay

1 U of yeast topoisomerase II will relax 0.5 μg of supercoiled pBR322 when incubated in Assay Buffer in a total reaction volume of 30 μl at 30 °C for 30 minutes. Gels are run in the absence of ethidium bromide or chloroquine.



## Quality Control

- 1) Purity: Yeast Topoisomerase II is purified to > 95 % purity as judged by SDS-polyacrylamide gel electrophoresis.
- 2) Tests for Yeast Topoisomerase I contamination by looking for relaxation of sc pBR322 under topoisomerase I assay conditions were negative.